



SPECS

- PN: HE-1256 and HE-1256/2
- Rebuilt to Zero-Hour Specs
- Compatible for Aquilion 4-Slice thru PRIME (Gen. 1) CT platforms

OUR HEAT EXCHANGERS VS. USED

- Completely Rebuilt and Tested
 Here
- 90-Day Warranty
- Expert Technical Support
- Ask about our convenient *Heat Exchanger Installation Service*

Call for details today

INNOVATION, CONVENIENCE & SAVINGS FROM THE LEADER IN TOSHIBA CT SERVICE SUPPORT.

As a leader in Toshiba Aquilion CT parts and training support, we are proud to offer Heat Exchangers that are completely rebuilt to zerohour specs. The process for heat exchangers at Richardson Electronics Healthcare is very similar to that of the OEM process at less than one fourth the OEM cost.

- All heat exchangers are completely disassembled
- Hardware replaced: Screws, washers, hoses, the quick disconnect fitting, and accumulator diaphragm
- The radiator, pump, and accumulator are chemically flushed to clean the interior surfaces
- The pressure relief valve is rebuilt.

LEARN MORE

This innovation is an example of how we draw upon our unique experience and engineering capabilities to support your CT service efforts. We have the industry's strongest team dedicated to listening to—and solving—your challenges.

- The pump, radiator, fans, pressure switch, and thermal switch are tested (If they do not meet specification they are replaced)
- Cabinet is repainted, if necessary
- The unit is reassembled, filled and tested.
- In addition, we will soon be able to apply our very own UL Safety Mark

Visit our website to learn more about our team, QA3-tested parts inventory, P3 Program, and the best Toshiba CT service training available.

LEARN MORE www.rellhealthcare.com