




RECYCLE PASSPORT

X-RAY TUBE HOUSING ASSEMBLIES

Products:	X-Ray Tube Housing Assembly
Models:	Alta750G
Manufacturer:	Richardson Electronics 40W267 Keslinger Road LaFox, Illinois 60147 USA
Electronic Information:	www.rellhealthcare.com
Total Weight:	60.5 kg.

Recycle Information

Special attention	Items	Location
	Product is heavy in nature.	
	Product contains high-vacuum sub-assembly.	
	Product is intended to return to the manufacturer for proper disposal.	
	Contact Richardson Electronics for proper handling and disposal.	
Fluids 	Items	Location
	Mineral Insulating Oil (No PCBs)	Housing
Hazardous 	Substances	Location
	Lead (for radiation shielding)	Housing
	Beryllium (Be)	X-Ray Tube Window Material

